



Grounding Pad Manufacturer	ConMed Brand Name	Skintact Order Number	ConMed Order Number	3M Order Number	Valleylab Brand Name	Valleylab Order Number
<b>ConMed</b>						
ConMed	Macrolyte	no Skintact plate, because single plate (no safety system)	400-2100	9135	Polyhesive II	E7506
ConMed	Macrolyte	no Skintact plate, because single plate (no safety system)	400-2349	9130		
ConMed	Macrolyte	no Skintact plate, because single plate (no safety system)	400-2700	9130F		
ConMed	Macrolyte	RS35A (also RO21 possible)	402-2204	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
ConMed	Macrolyte	RS35A (also RO21 possible)	402-2204-15	9165L	Polyhesive II	E7507-DB or E7510-25DE
ConMed	Macrolyte	RO21	402-2508	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
ConMed	Macrolyte	RS25	402-2750	9160 or 1180	Polyhesive II	E7509
ConMed	Macrolyte	RS25	402-2800	9160F or 1180	Polyhesive II	E7509B
ConMed	Macrolyte	no Skintact plate, because single plate (no safety system)	425-2200	9135		
ConMed	Macrolyte	no Skintact plate, because single plate (no safety system)	425-2725	9130		
ConMed	Macrolyte	no Skintact plate, because neonatal plate	440-2400	No 3M Equivalent		
ConMed	Macrolyte	no Skintact plate, because single plate (no safety system)	450-2300	No 3M Equivalent		
ConMed	SureFit	RO21	410-2000	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
ConMed	SureFit	RO21	410-2100	9165L	Polyhesive II	E7507-DB or E7510-25DE
ConMed	SureFit	RO21A30	410-2200	9160F or 1180	Polyhesive II	E7509B
ConMed	SureFit	RO21A30	410-2400	9160 or 1180	Polyhesive II	E7509
ConMed	ThermoGard	RO21 OR RS35A	51-7310	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
ConMed	ThermoGard	RO21A30 or RS25	51-7410	9160 or 1180	Polyhesive II	E7509
ConMed	ThermoGard	RO21	51-7710	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
ConMed	ThermoGard	no Skintact plate, because single plate (no safety system)	51-7810	9135		
ConMed	ThermoGard	no Skintact plate, because single plate (no safety system)	51-7910	9135		
ConMed	ThermoGard	RO21 OR RS35A	7-382	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
ConMed	ThermoGard	RO21	7-383	9165 or 1179/81	Polyhesive II	E7507 or E7510-25
<b>Valleylab</b>						
Valleylab	Polyhesive II	no Skintact plate, because neonatal plate		No 3M Equivalent		E7512