

# Stainless Steel Tables

## MCM520 Series

### Utility Tables

Utility Table with 3-sided guardrail and lower drawers. Standard size is 16"W x 20"L and 34"H, custom sizes available. 2 locking caster on front to secure unit in place. Locking top drawer available. Fully welded construction and made to OR specifications.

- > Drawers slide freely
- > No ripped gloves
- > Custom options available, get what you want


MCM520	One Drawer and Lower Shelf
MCM521	Two Drawers and Lower Shelf
MCM522	Three Drawer and Lower Shelf
MCM523	Four Drawers and Lower Shelf



If you want to make sure something doesn't fall off, we can add a **guardrail** to any table.



Need more storage? A drawer or **drawers** can be added to any table.



3" neoprene scuff resistant swivel are a standard feature. Locking **casters** available by request.